

# LUXSMART™

P R E L O A D E D

## PREMIUM HYDROPHOBIC IOL

For your daily  
range of vision



CATARACT



LASER



RETINA


**BAUSCH+LOMB**

See better. Live better.


## 70s ARE THE NEW 40s

The ESCRS Functional Vision Working Group reported that Europeans who are 55 years and older spend at least **6 hours per day on leisure activities**<sup>1</sup>, including playing games and computer use, relaxing/thinking, reading, watching television, socializing and communicating, participating in exercise, recreation, and other activities, including travel.


**Figure 1. People aged ≥ 50 years old spending at least 3 hours per week on physical activity outside work**



**Figure 2. People aged 65-74 years who have never used a computer, 2008 vs. 2017**



**Figure 3. Time spent on shopping and personal services, > 65 years old**



Source:



Besides leisure activities, several working distances are also needed for performing other common daily tasks, such as cooking, seeing the speedometer in a car, or walking on uneven ground.

**Figure 4. Average time that consumers who cook at home spend cooking each week (hours)**

	Germany	France	Italy	United Kingdom	Poland
Total	5.4	5.5	7.1	5.9	6.1
Women	6.5	6.7	8.8	6.3	7.6
Men	4.3	4.2	5.3	5.4	4.5
Aged 15-19	4.1	3.3	4.7	4.3	3.8
Aged 20-29	4.3	4.8	6.9	5.4	5.3
Aged 30-39	5.5	5.1	7.5	5.7	6.5
Aged 40-49	5.4	5.8	7.6	5.9	6.5
Aged 50-59	6.3	6.2	7.5	6.4	9.3
Aged 60 plus	6.4	6.7	7.0	6.5	6.9

Source:



## OPTICAL CONCEPT

### PURE REFRACTIVE OPTICS (PRO) Technology

With no diffractive optical profile; the IOL\* has a refractive surface across the entire optical diameter

#### PERIPHERY

Refractive aspheric surface

#### ELONGATED FOCUS CENTER

2 mm center with combination of 4<sup>th</sup> and 6<sup>th</sup> orders of spherical aberration of **opposite signs**

#### PATENTED TRANSITION ZONE

Transition zone designed to smoothly decrease the optic vergence from the center to the periphery

Transition designed to take part of the 4<sup>th</sup> and 6<sup>th</sup> orders of spherical aberration management

Transition designed to control the trajectory of light rays to ensure no light is outside the range of vision (no light loss)

\*IOL: Intraocular lens


## The Area under the Modulation Transfer Function (MTFa) and its relationship with the Visual Acuity

The MTFa is an objective in vitro MFT-based metric to assess the optical quality of an intraocular lens: the larger the MTFa value, the better the IOL optical quality

As opposed to MTF at single spatial frequency, the MTFa is the area under the MTF curve calculated from 0 to 50 cycles/mm.

Studies<sup>2,3,4</sup> have shown high correlation between MTFa and clinical visual acuity, so that it can be used to predict the visual performance at different levels of focus of pseudophakic patients.


**Figure 5. LuxSmart™ experimental Through-focus MTFa and predicted defocus range<sup>5</sup>**




For defocus values where the MTFa value is  $\geq 12$  (dotted line), the expected visual acuity would be around 0.2 log MAR (required for driving license)

**Figure 6. Depth of Focus comparison of experimental Through-focus MTFa and predicted defocus range for LuxSmart™ (Bausch + Lomb) and Acrysof™ IQ Vivity™ (Alcon)<sup>5</sup>**

For defocus values where the MTFa value is  $\geq 20$ , the expected visual acuity would be around 0.0 logMAR.



**Figure 7. Pinhole images and halos for LuxSmart™ (Bausch + Lomb) and Acrysof™ IQ Vivity™ (Alcon) at distance (top) and intermediate (+1.50 D) focus (bottom) at 4.5 mm pupil. Images are displayed in logarithmic scale for visualization purposes<sup>5</sup>**



Images of a pinhole object obtained at the distance focus of each lens with pupils of 4.5 mm. The images are displayed in logarithmic grayscale. The pinhole is a small but extended object which subtends an angle with respect to the model eye similar to the angle subtended by a car headlight of 10 cm observed at 100 m.

A double halo structure has an inner part with higher intensity due to the overlapping of the intermediate and distance defocused contributions

<sup>2</sup>. Visual acuity of pseudophakic patients predicted from in-vitro measurements of intraocular lenses with different design. F. Vega et al. Biomed. Opt. Express 9, 4893-4906 (2018).


<sup>3</sup>. Preclinical metrics to predict through-focus visual acuity for pseudophakic patients. A. Alarcon et al. Biomed. Opt. Express 7, 1877-1888 (2016).

<sup>4</sup>. Equivalence of two optical quality metrics to predict the visual acuity of multifocal pseudophakic patients. J. Armengol et al. Biomed. Opt. Express 11, 2818-2829 (2020)

<sup>5</sup>. Comparative optical bench analysis of a new extended range of vision intraocular lens. Juan Antonio Azor, Fidel Vega, Jesus Armengol, Maria S. Millan Grupo de Optica Aplicada y Procesado de Imagen (GOAPI). Department of Optics and Optometry Universitat Politecnica de Catalunya BARCELONATECH


## FOR OPTIMIZED EFFECTIVENESS AGAINST PCO\*

LuxSmart™ has a 360° continuous square edge on the posterior surface **to reduce incidence of posterior capsule opacification** in preventing epithelial lens cell migration under the IOL optic.<sup>6</sup>



Nixon and Woodcock<sup>7</sup> demonstrated that a **continuous 360° square edge had significantly less PCO than a square edge that was interrupted at the optic-haptic junction.**

## PROTECTION FROM UV LIGHT



**Figure 8.** Spectral transmission curves of LuxSmart™ and LuxSmart™ Crystal. The continuous vertical line marks the separation (380 nm) between the ultraviolet band and the visible spectrum.

\*PCO: Posterior capsule opacification

6. BAUSCH + LOMB data on file: RD-R-015. Measurement of sharp edge.

7. Nixon DR, Woodcock MG. Pattern of posterior capsule opacification models 2 years postoperatively with 2 single-piece acrylic intraocular lenses. J Cataract Refract Surg 2010; 36:929-934

## PLATFORM STABILITY

The shape of the LuxSmart™ has been designed to optimize its post-operative behavior in the capsular bag.

IOLs with a similar 4-point fixation haptic design have shown:

- › To have **good centration**<sup>8</sup>
- › To have similar **postoperative performances in terms of CDVA, inflammation and PCO** compared with the C-loop design<sup>8</sup>
- › To have **rotational stability**. 90 % of lenses rotate less than 5 degrees at 6 months<sup>9</sup>
- › To be **stable in the eye** and even suitable for the application of a toric surface to correct corneal astigmatism<sup>10</sup>

Orientation features of the LuxSmart™ IOL have been designed close to the optic edge **to facilitate visualization, specially in case of constricted iris.**

8. Mingels, A., Koch, J., Lommatzsch, A. et al. Comparison of two acrylic intraocular lenses with different haptic designs in patients with combined phacoemulsification and pars plana vitrectomy. Eye 21:1379-1383 (2007).


9. Kwartz J, Edwards K Evaluation of the long-term rotational stability of single-piece, acrylic intraocular lenses. British Journal of Ophthalmology 2010;94:1003-1006

10. Buckhurst, Phillip J., Wolffsohn, James S. PhD; Naroo, Shehzad A. PhD; Davies, Leon N. PhD Rotational and centration stability of an aspheric intraocular lens with a simulated toric design, Journal of Cataract & Refractive Surgery, September 2010 - Volume 36 - Issue 9 - p 1523-1528

## SINGLE STEP FULLY PRELOADED INJECTION

LuxSmart™ and LuxSmart™ Crystal are only available in a preloaded version, taking the advantage of:

- › **Less risk of IOL damage and mishandling<sup>11</sup>**
- › Usage of preloaded injection systems, which have been shown to produce **faster and more predictable IOL delivery<sup>11</sup>** with less wound stretching<sup>12</sup>
- › It is thought that during the next several years, use of **preloaded** disposable injectors is expected to **grow** and may well represent the **industry's future<sup>12,13</sup>**



**LUXSMART™**



**LUXSMART™ CRYSTAL**

## TECHNICAL SPECIFICATIONS

### MATERIAL

<b>Material:</b>	Acrylic hydrophobic
<b>Overall diameter:</b>	11.00 mm
<b>Optic diameter:</b>	6.00 mm
<b>Platform design:</b>	Single piece, 4 fixation points and 360° posterior square-edges
<b>Optical design:</b>	Asphericity modulation design with the combination of 4 <sup>th</sup> and 6 <sup>th</sup> orders of spherical aberration of <b>opposite signs</b>
<b>Haptics angulation:</b>	0°
<b>Light Filter:</b>	<b>LuxSmart™ Crystal:</b> UV filter <b>LuxSmart™:</b> UV and violet filters
<b>Dioptric range:</b>	From 0.00 D to +10.00 D (1.00 D steps) From +10.00 D to +34.00 D (0.50 D steps)
<b>Refractive index:</b>	1.54 at 35°
<b>Orientation features:</b>	Top right and bottom left

### DELIVERY SYSTEM


**Fully preloaded system with push injection:** Accuject™ Pro  
**Recommended incision size:** ≥ 2.2 mm (wound assisted technique)



### CONSTANTS\*

**OPTICAL CONSTANTS**  
**SRK/T Constant A:** 118.5  
**ACD:** 5.23  
**Surgeon factor:** 1.48

**Haigis:** a<sub>0</sub>: 1.045 / a<sub>1</sub>: 0.4 / a<sub>2</sub>: 0.1  
**EVO Constant A:** 118.5  
**Barrett:** Constant A: 118.4 / Lens Factor: 1.57  
**Hill-RBF 2.0 Constant A:** 118.32



Scan the code to access a real implantation video.  
*Courtesy of Dr. Hoerster, Germany*

<sup>11</sup> Chung B, Lee H, Choi M, Seo KY, Kim EK, Kim TI. Preloaded and non-preloaded intraocular lens delivery system and characteristics: human and porcine eyes trial. Int J Ophthalmol 2018;11(1):6-11  
<sup>12</sup> Mencucci R, Favuzza E, Salvatici MC, Spadea L, Allen D. Corneal incision architecture after IOL implantation with three different injectors: an environmental scanning electron microscopy study. Int Ophthalmol. Published online: 01 February 2018. <https://doi.org/10.1007/s10792-018-0825-2>

\*Constants are estimates only. It is recommended that each surgeon develops their own values.

# LUXSMART™

P R E L O A D E D



@BauschSurgical



Bausch + Lomb Surgical

[www.bauschsurgical.eu](http://www.bauschsurgical.eu)

CE  
0197

© 2021 Bausch + Lomb Incorporated.  
All rights reserved. ®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.  
All other brand/product names are trademarks of the respective owners.  
For healthcare professionals only, please refer to the instructions for use.

EMEA\_SU\_B\_LuxS\_21\_002



CATARACT



LASER



RETINA

**BAUSCH + LOMB**  
See better. Live better.