

BAUSCH+LOMB

Stellaris[®] *PC*

Vision Enhancement System

Product Catalogue
Fully integrated
Procedural versatility
Even more intuitive

choosestellaris.com



Stellaris[®] PC

- Complete combined surgical system
- Enhanced procedural versatility
- Enhanced interface for improved ease-of-use

Table of Contents

4	Hardware and Accessories
6	Combined Procedural Packs
8	Posterior Procedural Packs
10	Standalone Items
11	VitreoRetinal Tamponades
12	Laser Disposable and Accessories
14	Single Use VitreoRetinal Instruments
16	Anterior Phaco Packs
18	Phaco Needles & Sleeves
19	Cataract Single Use instruments



BL1433

Stellaris PC Combined System
with 1 Xenon + 1 Xenon-Mercury
Lamp Assembly

BL14304

Stellaris PC Next Generation
Combined System with an integrated
laser with 1 Xenon + 1 Xenon Mercury
Lamp Assembly

BL14334

Stellaris PC Next Generation
Combined System (without the
integrate laser) with 1 Xenon + 1
Xenon-Mercury Lamp Assembly



BL2295 Integrated Foot Pedal

BL4390 Foot Control Battery

BL4391 Foot Control Wall Charger

BL4393 Foot Control Charging Cradle

BL4394 Foot Control Power Cable

BL4363 Stellaris PC Zero Level Bottle Hanger -
Dual Bottle



BL3270 Fragmentation Handpiece



= 1



BL3170 Stellaris Phaco Handpiece



= 1



BL3320 Xenon Lamp Replacement



= 1

BL3321 Mercury Lamp Replacement



= 1

BL5275 20GA Fragmentation Needle 15 deg
single use

= 6



BL2272 Laser Upgrade kit



= 1

BL3255EN Stellaris PC Manual V4 English



= 1



BL4364 Stellaris PC/Stellaris Mayo Tray Adapter



= 1



Please contact your local representative to order a Mayo Tray.

06 / Stellaris PC Product Catalogue

Combined Procedural Packs with AFI (Air Forced Infusion)

BL5520W

20G Combined Pack

Pack Includes:

- A 20G Vitrectomy Cutter
- B Scleral plug package - 19 & 20
- C 20G Wide angle light pipe
- D 20G MVR knife



BL5523WV

23G Combined Pack

Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



BL5525WV

25G Combined Pack


Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Irrigation sleeves x 2
- 11 Needle wrench
- 12 Pre-piercing spike assembly



2

3

4



5

6

7



8

9

10

11

12

BL5320W**20G Combined Pack**

Pack Includes:

- A 20G Vitrectomy Cutter
- B Scleral plug package - 19 & 20
- C 20G Wide angle light pipe
- D 20G MVR knife



STERILE



B



C

D

BL5323WV**23G Combined Pack**

Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



STERILE



B



C



D

BL5325WV**25G Combined Pack**

Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



STERILE



B




C



D

All the above packs come with the items below.


- 1 Cassette with combined gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Irrigation sleeves x 2
- 11 Needle wrench



5

6

7



2

3



4



8



9



10

11



08 / Stellaris PC Product Catalogue

Posterior Procedural Packs with AFI (Air Forced Infusion)

BL5420W

20G Posterior Pack

Pack Includes:

- A 20G Vitrectomy Cutter
- B Scleral plug package - 19 & 20
- C 20G Wide angle light pipe
- D 20G MVR knife



BL5423WV

23G Posterior Pack

Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



BL5425WV

25G Posterior Pack


Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Pre-piercing spike assembly



2

3

4



5

6

7



8

9

10

BL5220W

20G Posterior Pack

Pack Includes:

- A 20G Vitrectomy Cutter
- B Scleral plug package - 19 & 20
- C 20G Wide angle light pipe
- D 20G MVR knife



C

D

BL5223WV

23G Posterior Pack

Pack Includes:

- A 23G Vitrectomy Cutter
- B 23G Wide angle light pipe
- C 23G Valved ESA System
- D 23G Scleral plugs



B



C

D

BL5225WV

25G Posterior Pack

Pack Includes:

- A 25G Vitrectomy Cutter
- B 25G Wide angle light pipe
- C 25G Valved ESA System
- D 25G Scleral plugs



B




C

D

All the above packs come with the items below.

- 1 Cassette with gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler




1

5

6

7



2

3



4



8



9



Focal

BL5720 20G Focal Light Pipe



BL5723 23G Focal Light Pipe



BL5725 25G Focal Light Pipe



Wide Angle

BL5820 20G Wide Angle Light Pipe



BL5823 23G Wide Angle Light Pipe



BL5825 25G Wide Angle Light Pipe



Cutters

BL5620 20G VIT Cutter



BL5623 23G VIT Cutter




BL5625 25G VIT Cutter



Valved trocars

BL5281 23G ESA kit-valved (kit of 3)



BL5282 25G ESA kit-valved (kit of 3)



BL5283 23G Single Trocar-Valved



BL5284 23G Replacement Valves



BL5285 25G Single Trocar-Valved



BL5286 25G Replacement Valves



Standard trocars

CX9123 23G ESA SYSTEM (single ESA)



CX9626 25G ESA SYSTEM (kit of 3)



New ESA valved cannula system
delivers improved performance and
enhanced surgeon control



VRL100 DK-Line® 5ml

STERILE (X)



= 1

VRL110 DK-Line® 7ml

STERILE (X)



= 1



VRL700 Oxane® HD, 10ml Syringe

STERILE (X)



= 1



VRL200 Okta-Line® 5ml

STERILE (X)



= 1



CX5710 Viscous Fluid Extraction Kits

STERILE (X)



= 6



VRL500 Oxane® 1300, 10ml Syringe

STERILE (X)



= 1

VRL600 Oxane® 5700, 10ml Syringe

STERILE (X)



= 1

CX5720 Viscous Fluid Injection Kits for
Oxane range

STERILE (X)



= 6



Stellaris® PC

With fully integrated 532nm laser

Introducing the next generation Stellaris PC. The ideal choice for cataract, VR and combined procedures, evolution driven by surgeon feedback.

Stellaris PC now features fully integrated laser. Upgrade compatible with existing Stellaris PC models.

BL2272 Laser Upgrade kit



= 1

Laser module
Laser mounting brackets
E-stop and cabling
Key switch and cabling (2) keys
Laser drawer replacement cover
Laser power input module
Anterior module with laser support
Laser to power input module cabling
Laser optical patch cord
CAN and power cabling
Upgrade software kit
External new labels and one pair of goggles

BL2295 Integrated new design footpedal



= 1



Standard Laser Probe

BL5290 20G STD. Laser PROBE Straight
BL5293 23G STD. Laser PROBE Straight
BL5295 25G STD. Laser PROBE Straight



= 5



Curved Illuminated Laser Probe

BL5290LC 20G ILL. Laser PROBE Curved
BL5293LC 23G ILL. Laser PROBE Curved
BL5295LC 25G ILL. Laser PROBE Curved



= 5



Curved Standard Laser Probe

BL5290C 20G STD. Laser PROBE Curved
BL5293C 23G STD. Laser PROBE Curved
BL5295C 25G STD. Laser PROBE Curved



= 5




Standard Aspirating Laser Probe

BL5290ASPH 20G Hard Tip asp. Laser PROBE
BL5293ASPH 23G Hard Tip asp. Laser PROBE



= 5



Illuminated Laser Probe

BL5290L 20G ILL. Laser PROBE Straight
BL5293L 23G ILL. Laser PROBE Straight
BL5295L 25G ILL. Laser PROBE Straight



= 5



Soft-tip Aspirating Laser Probe


BL5290ASP 20G Soft Tip asp. Laser PROBE S
BL5293ASP 23G Soft Tip asp. Laser PROBE S



= 5



Laser filters



BL3236 Zeiss Fixed ESF



BL3240 Zeiss 2 position filter ESF



BL3245 Moeller Weidel Fixed ESF



BL3241 Leica 2 position filter ESF



BL3247 Leica 4-way filter



BL3237 Leica Fixed ESF



BL3239 Laser filter ESF Mount Topcon



BL3242 Adaptor for fitting of 2 position filters



BL3246 LIO Replacement Bulb



BL3234 Laser On/Off Keys (1 set)



BL3235 Laser Safety Goggles



BL3244 Interlock bypass connector kit

- Microscope filter override key
- Microscope filter connector
- Laser room interlock connector



BL2298 Dual LIO



BL2296 Auxiliary脚踏 for LIO use



Endgripping Forceps

65023S-12 23G STERILE   = 12

65025S-12 25G STERILE   = 12



Asymmetric Peeling Forceps

SUI03-12 23G STERILE   = 12

SUI05-12 25G STERILE   = 12



ILM Forceps

SUI09-12 20G STERILE   = 12

SUI07-12 23G STERILE   = 12

SUI08-12 25G STERILE   = 12



Serrated Forceps

SUI10-12 20G STERILE   = 12

SUI11-12 23G STERILE   = 12

SUI12-12 25G STERILE   = 12



Membrane Picks

SUK01 23G STERILE   = 8



SUK02 25G STERILE   = 8



Bipolar Eraser

E7928 Scleral Bipolar Eraser, Lemos type connector STERILE   = 12



20G

E7929 20G Bipolar Eraser, Tapered, Lemos type connector STERILE   = 12



23G

E7930 23G Bipolar Eraser, Tapered, Lemos type connector STERILE   = 12



E7931 23G Bipolar Eraser, Tapered, 2-pin connector STERILE   = 12



25G

E7924 25G Bipolar Eraser, Tapered, 2-pin connector STERILE   = 12



E7925 25G Bipolar Eraser, Blunt, 2-pin connector STERILE   = 12



E7926 25G Bipolar Eraser, Tapered, Lemos type connector STERILE   = 12



E7927 25G Bipolar Eraser, Blunt, Lemos type connector STERILE   = 12



Aspiration Handpieces

SUJ03 23G Passive Aspiration Handpiece, 1.5mm soft tip cannula

STERILE (X)  = 12



SUJ02 25G Passive Aspiration Handpiece, 1.5mm soft tip cannula

STERILE (X)  = 12



SUJ06 23G Passive Aspiration Handpiece with blunt needle

STERILE (X)  = 12



SUJ07 25G Passive Aspiration Handpiece with blunt needle

STERILE (X)  = 12



SUJ08 20G Backflush/Extrusion Handpiece, blunt tip, grey color code ring

STERILE (X)  = 6



SUJ09 20G Backflush/Extrusion handpiece. Soft tip brush, grey color code ring. Not compatible with valved trocar/cannula system

STERILE (X)  = 6



SUJ10 20G Backflush/Extrusion handpiece. Soft tip, grey color code ring. Not compatible with valved trocar/cannula system

STERILE (X)  = 6



SUJ11 23G Backflush/Extrusion handpiece. Blunt tip, green color code ring.

STERILE (X)  = 6



SUJ12 23G Backflush/Extrusion handpiece. Soft tip brush, green color code ring. Not compatible with valved trocar/cannula system

STERILE (X)  = 6



SUJ13 23G Backflush/Extrusion handpiece. Soft tip, green color code ring. Not compatible with valved trocar/cannula system

STERILE (X)  = 6



SUJ14 25G Backflush/Extrusion handpiece. Blunt tip, blue color code ring.

STERILE (X)  = 6



MVR Knives

Laseredge®




E7519 19G MVR Blade
A: 31.8mm B: 1.5mm C: 3.1mm

STERILE (X)  = 6

E7520 20G MVR Blade
A: 31.8mm B: 1.2mm C: 2.4mm

E7523 23G MVR Blade
A: 31.8mm B: 0.6mm C: 1.6mm

E7525 25G MVR Blade
A: 31.8mm B: 0.5mm C: 1.0mm



BL5110


Stellaris Premium Vacuum
Phacoemulsification Pack



Pack Includes:

- A VFM cassette assembly and I/A set
- B Auxillary drape
- C Screen drape

- D Remote control drape
- E Irrigation sleeves x 2
- F Plug
- G Test chamber
- H Needle wrench



BL5111


Stellaris Basic Vacuum
Phacoemulsification Pack



Pack Includes:

- A VFM cassette assembly and I/A set
- B Irrigation sleeve
- C Plug

- D Test chamber
- E Needle wrench



BL5112


Stellaris Stable Chamber Vacuum Pack



Pack Includes:

- A VFM cassette assembly and I/A set
- B Auxillary drape
- C Screen drape
- D Remote control drape

- E Irrigation sleeves x 2
- F Plug
- G Test chamber
- H Needle wrench
- I StableChamber



BL5113


Stellaris StableChamber™ Phaco Pack
with MICS straight needle



Pack Includes:

- A VFM cassette assembly and I/A set
- B Auxillary drape
- C Screen drape
- D Remote control drape

- E Irrigation sleeves x 2
- F Plug
- G Test chamber
- H Needle wrench
- I StableChamber
- J MICS straight needle



BL5113A


Stellaris StableChamber™ Phaco Pack
with MICS angled needle



Pack Includes:

- A VFM cassette assembly and I/A set
- B Auxillary drape
- C Screen drape
- D Remote control drape


- E Irrigation sleeves x 2
- F Plug
- G Test chamber
- H Needle wrench
- I StableChamber
- J MICS angled needle



BL5114
 Stellaris Premium Vacuum Phaco Pack
 with MICS straight needle




- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I MICS straight needle



BL5114A
 Stellaris Premium Vacuum Phaco Pack
 with MICS angled needle



- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I MICS angled needle



BL5612
 20GA High Speed
 (30-5000cpm) Anterior
 Vitrectomy Cutter Pack



BL5280 Dual Infusion Kit



D4600A DigiFlow™ Tubing



Phaco tips for MICS® Surgery

RIS: Recommended Incision Size in Coaxial

OD: Outer Diameter

ID: Inner Diameter

Single Use

BL3318S MICS® needle 1,8

STERILE   = 6

BL3318AS MICS® needle 1,8 Angled 30°

STERILE   = 6

RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm

Reusable limited re-use (10 cases)

BL3318 MICS® needle 1,8

 = 6

BL3318A MICS® needle 1,8 Angled 30°

 = 6

RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm

BL3118 Kit of 2 SU MICS Sleeves

STERILE   = 12

For Standard Surgery

RIS: Recommended Incision Size in Coaxial

OD: Outer Diameter

ID: Inner Diameter

Single Use

DP8130S Standard needle 30° 19G

STERILE   = 6

DP8145S Standard needle 45° 19G

STERILE   = 6

RIS: 3.00-3.20mm OD: 1.70mm ID: 0.90mm

DP8230S Microflow needle 30°

STERILE   = 6

DP8245S Microflow needle 45°

STERILE   = 6

RIS: 2.50-2.65mm OD: 1.07mm ID: 0.90-0.51mm

DP8345S Microflow Plus needle 45°

STERILE   = 6

DP8330S Microflow Plus needle 30°

STERILE   = 6

RIS: 2.75-2.85mm OD: 1.17mm ID: 0.98-0.61mm

DP8430S Thin Tip needle 30° 20G

STERILE   = 6

RIS: 2.5-2.65mm OD: 0.90mm ID: 0.67mm

DP8730S Standard needle 10°

STERILE   = 6

RIS: 2.2mm OD: 0.90mm ID: 0.70-0.50mm

Reusable limited re-use (10 cases)

DP8130 Standard needle 30° 19G

 = 6

DP8145 Standard needle 45° 19G

 = 6

RIS: 3.00-3.20mm OD: 1.07mm ID: 0.90mm

DP8200 Microflow needle 0°

 = 6

DP8215 Microflow needle 15°

 = 6

DP8230 Microflow needle 30°

 = 6

DP8230A Microflow needle 30° Angled

 = 6

DP8245 Microflow needle 45°

 = 6

RIS: 2.50-2.65mm OD: 1.07mm ID: 0.90-0.51mm

DP8300 Microflow Plus needle 0°

 = 6

DP8315 Microflow Plus needle 15°

 = 6

DP8330 Microflow Plus needle 30°

 = 6DP8330A Microflow Plus needle 30°
Angled = 6

DP8345 Microflow Plus needle 45°

 = 6

RIS: 2.75-2.85mm OD: 1.17mm ID: 0.98-0.61mm

DP8415 Thin Tip needle 15° 20G

 = 6

DP8430 Thin Tip needle 30° 20G

 = 6

RIS: 2.5-2.65mm OD: 0.90mm ID: 0.67mm

DP5531 kit of 2 SU Sleeves 2,2 surgery

STERILE   = 12

DP5501 Standard kit of 2 SU Sleeves

STERILE   = 12

Manipulators

SUF01 Koch Stop & Chop Manipulator STERILE   = 12



SUF02 Sinskey Hook STERILE   = 12



SUF03 Phaco Spatula STERILE   = 12



SUF04 Drysdale Nucleus Manipulator STERILE   = 12



SUF05 Horizontal Manipulator with Y-Hook STERILE   = 12



SUF06 Vertical Manipulator with Y-Hook STERILE   = 12



Cannulas

82020S Corneal Irrigator STERILE   = 60



Speculum

SUH01 Barraquer Eye Speculum STERILE   = 24
available as Per Procedure Tray component part



SUH04 Eye Speculum, adjustable STERILE   = 6



Forceps

SUE01-12 Capsulorhexis Forceps STERILE   = 12



SUE02-12 Bonn Forceps STERILE   = 12



SUE03-12 Tying Forceps STERILE   = 12



SUE04-12 McPherson Tying Forceps, angled STERILE   = 12



SUE05-12 Pierse Tying Forceps STERILE   = 12



SUE06-12 MICST[™] Capsulorhexis Forceps, fine angled tips, angled shafts for 1.8mm incision. STERILE   = 12



SUE07-12 MICST[™] Capsulorhexis with markings STERILE   = 12




SUE08-12 MICST[™] Capsulorhexis with curved jaws STERILE   = 12



Laseredge®

Guarded



Clear Corneal Knives

Laseredge®

E7594 1.3mm Angled Double Bevel
A: 12.6mm B: 1.3mm C: 1.4mm

E7595 1.6mm Angled Double Bevel
A: 12.6mm B: 1.6mm C: 1.7mm

E7527 2.75mm Angled Double Bevel
A: 14.4mm B: 2.7mm C: 3.1mm

E7528 2.85mm Angled Double Bevel
A: 14.4mm B: 2.8mm C: 3.1mm

E7587A 2.9mm Angled Double Bevel
A: 14.4mm B: 2.9mm C: 3.2mm


E7529 3.0mm Angled Double Bevel
A: 14.4mm B: 3.0mm C: 3.2mm

E7532 3.2mm Angled Double Bevel
A: 14.4mm B: 3.2mm C: 3.6mm

Guarded

E7528G 2.85mm Angled Double Bevel
A: 18.5mm B: 2.8mm C: 3.1mm

STERILE   = 6



Slit Knives continued

Laseredge®

E7551A 2.65mm Angled Bevel Up
A: 14.6mm B: 2.6mm C: 3.9mm

E7551AT 2.65mm Angled Bevel Up
Tip Only A: 14.6mm B: 2.6mm C: 3.9mm

E7551ADB 2.65mm Angled Double Bevel
A: 14.6mm B: 2.6mm C: 3.9mm

E7559A 2.75mm Angled Bevel Up
A: 14.4mm B: 2.7mm C: 3.2mm

E7552ADB 2.80mm Angled Double Bevel
A: 14.4mm B: 2.8mm C: 3.4mm

E7553A 2.85mm Angled Bevel Up
A: 14.1mm B: 2.8mm C: 3.3mm

E7556A 3.0mm Angled Bevel UP
A: 14.4mm B: 3.0mm C: 3.9mm

E7556ADB 3.0mm Angled Double Bevel
A: 14.4mm B: 3.0mm C: 3.7mm

E7557A 3.2mm Angled Bevel UP
A: 14.4mm B: 3.2mm C: 3.9mm


E7557 3.2mm Straight
A: 14.4mm B: 3.2mm C: 3.9mm

Guarded

E7550AG 2.5mm Angled Bevel UP
A: 18.5mm B: 2.5mm C: 3.1mm

E7551AG 2.65mm Angled Bevel UP
A: 18.5mm B: 2.6mm C: 3.8mm

STERILE   = 6



Slit Knives

Laseredge®

E7598 1.4mm Angled Bevel Up
A: 13.6mm B: 1.4mm C: 1.4mm

E7544A 1.8mm Angled Bevel Up
A: 14.4mm B: 1.8mm C: 2.7mm

E7544ADB 1.8mm Angled Double Bevel
A: 14.4mm B: 1.8mm C: 2.2mm


E7548ADB 2.2mm Angled Double Bevel
A: 14.4mm B: 2.2mm C: 2.7mm

E7549ADB 2.4mm Angled Double Bevel
A: 14.4mm B: 2.4mm C: 2.9mm

E7550A 2.5mm Angled Bevel Up
A: 14.4mm B: 2.5mm C: 3.8mm

E7550ADB 2.6mm Angled Double Bevel
A: 14.4mm B: 2.6mm C: 3.3mm

STERILE   = 6




Stab Knives

Laseredge®

E7515 15°
A: 11.2mm B: 1.6mm C: 5.9mm

STERILE   = 6



Guarded


E7515G 15°
A: 14.0mm B: 1.6mm C: 5.9mm

Trapezoidal Knives

Laseredge®

E7599 1.5-1.7mm Angled
A: 14.4mm B: 1.7mm C: 2.0mm D: 1.5 E: 2.9

STERILE   = 6



E7596 1.5-2.0mm Angled
A: 14.4mm B: 2.0mm C: 1.8mm D: 1.5 E: 3.4

E7600 1.6-1.8mm Angled
A: 14.4mm B: 1.8mm C: 1.8mm D: 1.6 E: 2.9

E7602 1.8-2.2-2.8mm Angled
A: 16.8mm B: 2.2mm C: 1.5mm D: 1.8 E: 3.4

E7603 1.85-2.2mm Angled
A: 14.4mm B: 2.2mm C: 1.8mm D: 2.0 E: 3.7

E7601 2.0-2.2mm Angled
A: 14.4mm B: 2.2mm C: 2.3mm D: 2.0 E: 3.1

E7572 2.75-3.2mm Angled Bevel Up
A: 16.8mm B: 3.2mm C: 3.2mm D: 2.7 E: 4.6mm

E7582 2.85-3.2mm Angled Bevel Up
A: 16.8mm B: 3.2mm C: 3.4mm D: 2.8 E: 4.3

Guarded

E7601G 2.0-2.2mm Angled
A: 18.5mm B: 2.2mm C: 2.3mm D: 2.0 E: 3.1

Stellaris® System Coaxial I/A Handpieces

85794ST MICS® I/A handpiece curved
shaft, with sleeve for 1.8mm incision

STERILE   = 12



85795ST MICS® I/A handpiece 45°
angled tip, with sleeve for 1.8mm incision

STERILE   = 12



85783ST 17G (Ø1.5mm) I/A Handpiece,
curved shaft

STERILE   = 12



85784ST 17G (Ø1.5mm) I/A Handpiece
straight shaft 45° angled tip 0.3mm
aspiration port

STERILE   = 12



85786ST 16G (Ø1.65mm) I/A
Handpiece straight shaft, 90° angled tip,
0.3mm aspiration port

STERILE   = 12



Stellaris® System Capsule Guard I/A® Coaxial Handpieces

85910ST CapsuleGuard® I/A Handpiece
45° Tip - Incisions: 2.2mm-2.8mm

STERILE   = 12



85912ST CapsuleGuard® I/A Handpiece
Straight Tip - Incisions: 2.2mm-2.8mm

STERILE   = 12



85913ST CapsuleGuard® I/A Handpiece
Curved Tip - Incisions: 2.2mm-2.8mm

STERILE   = 12



85914ST CapsuleGuard® I/A Handpiece MICS® Straight Tip – Incisions: sub 20mm

STERILE



= 12



85915ST CapsuleGuard® I/A Handpiece MICS® 45° Tip – Incisions: sub 20mm

STERILE



= 12



Stellaris® System Biaxial I/A Handpieces Aspiration

85901S19 Aspiration handpiece 19G (Ø 1.07mm)

STERILE



= 12

85901S20 Aspiration handpiece 20G (Ø 0.89mm)

STERILE



= 12

85780S Aspiration handpiece 21G (Ø 0.8mm)

STERILE



= 12

85790S Aspiration handpiece 22G (Ø 0.7mm)

STERILE



= 12

85792S Aspiration handpiece 23G (Ø 0.6mm)

STERILE



= 12



curved shaft roughened tip

Stellaris® System Biaxial I/A Handpieces Irrigation

85902ST19 1.8mm B-MICS 19G (Ø 1.07mm) Irrigation handpiece

STERILE



= 12

85787ST 21G Irrigation handpiece (Ø 0.8mm)

STERILE



= 12

85806ST 22G Irrigation handpiece (Ø 0.7mm)

STERILE



= 12

85807ST 23G Irrigation handpiece (Ø 0.6mm)

STERILE



= 12



curved shaft smooth tip

Stellaris® System Irrigating Chopper for biaxial MICS

85903ST19 19G (Ø 1.07mm) Nagahara type tip, end out

STERILE



= 12



Stellaris® System Coaxial Zero Phaco Handpieces

SUP03ST Stellaris® System Coaxial Zero Phaco Handpiece

STERILE



= 12

Designed for lens removal following Femtosecond laser fragmentation without the use of ultrasonic energy. 30° beveled aspiration tube, pre-assembled with infusion sleeve. Recommended for standard incision >2.4mm.



SUP04ST Stellaris® System Coaxial Zero Phaco Handpiece 2.2mm

STERILE



= 12

Designed for lens removal following Femtosecond laser fragmentation without the use of ultrasonic energy. 30° beveled aspiration tube, pre-assembled with infusion sleeve. Recommended for 2.2mm incision.



Associated Marketing tools

Stellaris® PC: Next Generation
choosstellaris.com



Stellaris® PC Brochure



ESA Brochure

Stellaris® Brochure



Laser Probes flash card

Stellaris® PC: Next Generation

700+ Stellaris PC's installed in Europe

Try the Stellaris PC for yourself at our European Centres of Excellence




Why I choose Stellaris® PC:
“We asked for a footpedal with integrated laser controls, we have it. Dual linear operation, we have it. We wanted illumination in the anterior segment in cases of corneal opacity, it's done. You have listened and Stellaris has evolved”.

Dr Dov Weinberger
Head of Ophthalmology
Rabin Medical Centre
[Watch my full interview at
choosestellaris.com/weinberger](http://www.choosestellaris.com/weinberger)

For more information visit
www.choosestellaris.com

For more information on instruments visit
www.storzeye.eu
or email storz-instruments@bausch.com



Stellaris PC Vision Enhancement System. ®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. © Bausch & Lomb Incorporated.*

EMEA_SPCNG_CAT_2015

BAUSCH + LOMB
See better. Live better.