

Stellaris Elite

Vision Enhancement System

Stellaris ELITE™ Catalogue




BAUSCH+LOMB

Table of Contents

- 4 Vitrectomy Solutions for Any Practice
(Hardware and Accessories)
- 9 Illumination
- 10 Laser
- 12 Instrumentation
- 19 Vitreoretinal Tamponades
- 20 Anterior Capabilities
- 23 Product References





BL15455
Stellaris ELITE™
Combined system with integrated
laser with 1 Xenon + 1 Xenon Mercury
lamp assembly


BL14455
Stellaris ELITE™
Combined system (without the
integrated laser) with 1 Xenon + 1
Xenon-Mercury lamp assembly

BL2394 Stellaris ELITE™ Integrated Foot Pedal
BL4390 Foot Control Battery
BL4391 Foot Control Wall Charger
BL4393 Foot Control Charging Cradle
BL4394 Foot Control Power Cable



BL2300HD DMS (Digital Media System)  = 1



BL4363 Stellaris ELITE™ Zero Level Bottle
Hanger - Dual Bottle  = 1



BL3270 Fragmentation Handpiece



BL3170 Stellaris Phaco Handpiece



BL3320 Xenon Lamp Replacement



BL3321 Mercury Lamp Replacement



BL5275 20G Fragmentation Needle 15°
Single Use



BL5277 23G Fragmentation Needle 30°
Single Use



BL4364 Stellaris ELITE™/Stellaris Mayo Tray Adapter



Please contact your local representative to order a Mayo Tray.

NEW!

SE5523MVB / SE5523WVB

23G Combined Pack

Pack Includes:

- A 23G Bi-Blade® Vitrectomy cutter **•bi•Blade®**
- B 23G Wide angle or midfield light pipe
- C 23G Valved ESA system
- D 23G Scleral plugs



SE5525MVB / SE5525WVB

25G Combined Pack

Pack Includes:

- A 25G Bi-Blade® Vitrectomy cutter **•bi•Blade®**
- B 25G Wide angle or midfield light pipe
- C 25G Valved ESA system
- D 25G Scleral plugs



SE5527MVB

27G Combined Pack

Pack Includes:

- A 27G Bi-Blade® Vitrectomy cutter **•bi•Blade®**
- B 27G Midfield light pipe
- C 27G Valved ESA system





BL5280 **Dual Infusion Kit**



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second aspiration line
- 3 Infusion line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Test chamber
- 7 Male luer coupler
- 8 Irrigation MICS sleeves x 2
- 9 Priming cup
- 10 Needle wrench
- 11 Pre-piercing spike assembly
- 12 Second stopcock



NEW!

SE5423MVB / SE5423WVB

23G Posterior Pack

Pack Includes:

- A 23G Bi-Blade® Vitrectomy cutter **bi-Blade**
- B 23G Wide angle or midfield light pipe
- C 23G Valved ESA system
- D 23G Scleral plugs



SE5425MVB / SE5425WVB

25G Posterior Pack

Pack Includes:

- A 25G Bi-Blade® Vitrectomy cutter **bi-Blade**
- B 25G Wide angle or midfield light pipe
- C 25G Valved ESA system
- D 25G Scleral plugs



SE5427MVB

27G Posterior Pack



Pack Includes:

- A 27G Bi-Blade® Vitrectomy cutter **bi-Blade**
- B 27G Midfield light pipe
- C 27G Valved ESA system



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second aspiration line
- 3 Infusion line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming cup
- 7 Male luer coupler
- 8 Pre-piercing spike assembly
- 9 Second stopcock



Cutter for Stellaris Elite™ not compatible with Stellaris® PC

SE5623B 23G Bi-Blade® Vitrectomy cutter *

bi-Blade®



SE5625B 25G Bi-Blade® Vitrectomy cutter

bi-Blade®



SE5627B 27G Bi-Blade® Vitrectomy cutter

bi-Blade®



BL3394 Handle extension



Valved Trocars

BL5281 23G ESA Kit-valved (kit of 3)



BL5282 25G ESA Kit-valved (kit of 3)



BL5240 27G ESA Kit-valved (kit of 3)



BL5283 23G Single Trocar-Valved



BL5284 23G Replacement Valves



BL5285 25G Single Trocar-Valved



BL5286 25G Replacement Valves



BL5241 27G Single Trocar-Valved



Mid-Field

BL5823M 23G



 = 6
 STERILE 

BL5825M 25G



 = 6
 STERILE 

BL5827M 27G



 = 6
 STERILE 

Diffuse Wide Angle

BL 5823D 23G



 = 6
 STERILE 

BL 5825D 25G



 = 6
 STERILE 

BL 5827D 27G



 = 6
 STERILE 

Adjustable Chandelier

BL5234 23G

BL5235 25G



 = 5
 STERILE 

NEW!

25G Vivid Chandelier

BL5235S 25G



 = 6
 STERILE 

27G Vivid Chandelier *

BL5237S 27G



 = 6
 STERILE 

Stellaris ELITE™

With fully integrated laser

BL2272 Laser Upgrade Kit



- Laser module
- Laser mounting brackets
- E-stop and cabling
- Key switch and cabling (2) keys
- Laser drawer replacement cover
- Laser power input module
- Anterior module with laser support
- Laser to power input module cabling
- Laser optical patch cord
- CAN and power cabling
- Upgrade software kit
- External new labels Goggles

BL2394 Integrated Foot Pedal



Laser Filters



BL3236 Zeiss Fixed ESF



BL3240 Zeiss 2 Position Filter ESF



BL3245 Moeller Wedel Fixed ESF



BL3241 Leica 2 Position Filter ESF



BL3247 Leica 4-way Filter



BL3237 Leica Fixed ESF



BL3239 Laser Filter ESF Mount Topcon



BL3242 Adaptor for Fitting of 2 Position Filters



BL3246 LIO Replacement Bulb



BL3234 Laser On/Off Keys (1 set)



BL3235 Laser Safety Goggles



BL3244 Interlock Bypass Connector Kit

- Microscope filter override key
- Microscope filter connector
- Laser room interlock connector



BL2298 TruFocus LIO Premiere™ with LED illumination



Directional™ Laser Probe

- 55.26.20G
- 55.26.23 23G
- 55.26.25 25G
- 55.26.27 27G



Illuminated Laser Probes

Flexible Tapered Laser Probe

- 55.25.23 23G
- 55.25.25 25G



BL2296 Auxiliary Footpedal for LIO Use



Curved Illuminated Laser Probe



- BL5290LC 20G Illuminated Laser Probe Curved
- BL5293LC 23G Illuminated Laser Probe Curved
- BL5295LC 25G Illuminated Laser Probe Curved



Flexible Tapered Illuminated Laser Probe*

- 55.70.23P 23G
- 55.70.25P 25G



Illuminated Directional™ Laser Probe*

- 55.47.23P 23G
- 55.47.25P 25G



* Adaptor is required with Stellaris® Elite

Storz® Ophthalmic Instruments

- Fully disposable
- Ergonomic handle with actuation from two sides
- Non reflective instrument shaft
- Color coded handles for gauge identification
- 20G, 23G, 25G, 27G



Endgripping Forceps

65023S-12 23G
65025S-12 25G

STERILE   = 12

SUI13-12 27G



Asymmetric Peeling Forceps

SUI03-12 23G
SUI05-12 25G

STERILE   = 12



ILM Forceps

SUI09-12 20G
SUI07-12 23G
SUI08-12 25G

STERILE   = 12



Serrated Forceps

SUI10-12 20G
SUI11-12 23G
SUI12-12 25G
SUI14-12 27G

STERILE   = 12



Selected sterile single use instruments are
available as component of a **Per Procedure Tray**.

The Per Procedure Tray is a complete kit of sterile single use instruments selected by you. Assemble your individual kit of sterile single use instruments with the online Per Procedure Tray Configurator. For more information visit: www.storzeye.eu



Membrane Picks

SUK01 23G   = 8



SUK02 25G   = 8



CX9430 Reusable Bipolar Cord
with Lemoso Connector



Bipolar Erasers

E7928 Scleral Bipolar Eraser,
Lemoso Type Connector   = 12



20G

E7929 20G Bipolar Eraser, Tapered,
Lemoso Type Connector   = 12



23G

E7930 23G Bipolar Eraser, Tapered,
Lemoso Type Connector   = 12



25G

E7926 25G Bipolar Eraser, Tapered,
Lemoso Type Connector   = 12



27G

E7932 27G Bipolar Eraser, Tapered,
Lemoso Type Connector   = 12



Aspiration Handpieces

SUJ03 23G Passive Aspiration Handpiece. 1.5mm soft tip cannula. STERILE   = 12



SUJ02 25G Passive Aspiration Handpiece. 1.5mm soft tip cannula. STERILE   = 12



SUJ06 23G Passive Aspiration Handpiece. Blunt needle. STERILE   = 12





SUJ07 25G Passive Aspiration Handpiece. Blunt needle. STERILE   = 12





SUJ08 20G Backflush/Extrusion Handpiece. Blunt tip. Grey color code ring. STERILE   = 6



SUJ09 20G Backflush/Extrusion Handpiece. Soft tip brush. Grey color code ring. STERILE   = 6





SUJ10 20G Backflush/Extrusion Handpiece. Soft tip. Grey color code ring. STERILE   = 6



SUJ11 23G Backflush/Extrusion Handpiece. Blunt tip. Green color code ring. STERILE   = 6



SUJ12 23G Backflush/Extrusion Handpiece. Soft tip brush. Green color code ring. STERILE   = 6



SUJ13 23G Backflush/Extrusion Handpiece. Soft tip. Green color code ring. STERILE   = 6



SUJ14 25G Backflush/Extrusion Handpiece. Blunt tip. Blue color code ring. STERILE   = 6



SUJ17 27G Backflush/Extrusion Handpiece. Blunt tip. Magenta color code ring. STERILE   = 6



MVR Knives

LaseredgePlus®

LP7519 19G MVR Blade A: 1.37mm STERILE   = 6

LP7520 20G MVR Blade A: 1.2mm



Synergetics® Pinnacle 360°™

Pinnacle 360° is our high quality line of fully disposable, hand-held retinal surgical instruments.

- Offers wide variety of disposable tip options
- Fully disposable
- 360° actuation to approach the retina from any angle
- Color coded handles for gauge identification
- 20G, 23G, 25G, 27G, Foreign Body



Pinnacle 360°

20G



23G



25G



27G



Foreign Body



ILM Forceps

NEW!

Reddy End-Grasping Forceps

39.32.23PIN 23G
39.32.25PIN 25G

STERILE   = 5



Eckardt Forceps

39.08PIN 20G
39.08.23PIN 23G
39.08.25PIN 25G
39.08.27PIN 27G

STERILE   = 5



End-Grasping Forceps

09.12PIN 20G
09.12.23PIN 23G
09.12.25PIN 25G

STERILE   = 5



Fine Tip Eckardt Forceps

09.10.23PIN 23G
09.10.25PIN 25G

STERILE   = 5



Microserrated™ End-Grasping Forceps

39.12.23PIN 23G
39.12.25PIN 25G

STERILE   = 5



Tano Asymmetrical Micro Forceps

39.03.23PIN 23G
39.03.25PIN 25G
39.03.27PIN 27G

STERILE   = 5



Heavy Membrane Forceps

Pyramid Tip Forceps

11.04.23PIN 23G
11.04.25PIN 25G

STERILE   = 5



Serrated Jaw Forceps

11.06.23PIN 23G
11.06.25PIN 25G

STERILE   = 5



Dual Purpose Forceps

Apte Corset Forceps

39.07.23PIN 23G
39.07.25PIN 25G

STERILE   = 5



d™ Tewari Omnicep Forceps

39.15.23PIN 23G
39.15.25PIN 25G
39.15.27PIN 27G

STERILE   = 5



Foreign Body

BB Removal Forceps

11.40PIN 18G

STERILE   = 5



Scissors

Horizontal Scissors Curved

23.76PIN 20G
23.76.23PIN 23G
23.76.25PIN 25G

STERILE   = 5



Foreign Body Forceps

11.20PIN 19G

STERILE   = 5



Morris Vertical Scissors

10.24PIN 20G
10.24.23PIN 23G
10.24.25PIN 25G

STERILE   = 5



Foreign Body Magnet

12.21PIN 18G

STERILE   = 5



Myopic Forceps

Fine Tip Eckardt Forceps

09.10.23PIN-M 23G
09.10.25PIN-M 25G
• 37mm

STERILE   = 5



Pick Forceps

Pick Forceps

09.05.23PIN 23G
09.05.25PIN 25G

STERILE   = 5



Microserrated™ Tewari Omnicep Forceps

39.15.23PIN-M 23G
39.15.25PIN-M 25G
• 37mm

STERILE   = 5



Adeleberg Forceps

11.46.25PIN 25G

STERILE   = 5

DDMS™ Diamond Dusted Membrane Scrapers

Tano

20.04.20 20G
20.04.23 23G
20.04.25S (Stiff) 25G

STERILE   = 5

20.04.27 27G



Retractable Tano

20.16.20 20G
20.16.23 23G
20.16.25 25G

STERILE   = 3

STERILE   = 5

STERILE   = 5



Bipolar

Bipolar Cautery

53.00.23 23G
53.00.25 25G
53.00.27 27G

STERILE   = 6



Visco Dissector

Hubbard Visco Dissector

12.06.25 25G 28G Tip
• Injects viscoelastic for delaminating large membranes from the retina

STERILE   = 6



12.06.27 27G
• Injects viscoelastic for delaminating large membranes from the retina

STERILE   = 6



Backflush

D31.08.23
23G Silicon Tip active backflush

STERILE   = 5

D31.08.25
25G Silicon Tip active backflush

STERILE   = 5



Subretinal Cannula

12.03.25
Subretinal cannula 25G tip 41G

STERILE   = 5



Soft tip cannula

30.08.23
23G Soft tip cannula

STERILE   = 6

30.08.25
25G Soft tip cannula

STERILE   = 6

30.08.27
27G Soft tip cannula

STERILE   = 6



BL7600 Universal viscous fluid control pack

STERILE   = 6



VRL100 DK-Line® 5ml

STERILE   = 1

VRL110 DK-Line® 7ml

STERILE   = 1



BL7623 High Flow VFI Cannula 23G

STERILE   = 6



VRL200 Okta-Line® 5ml

STERILE   = 1

BL7625 High Flow VFI Cannula 25G

STERILE   = 6



VRL500 Oxane® 1300, 10ml Syringe

STERILE   = 1

VRL600 Oxane® 5700, 10ml Syringe

STERILE   = 1



VRL700 Oxane® HD, 10ml Syringe

STERILE   = 1



NEW!

BL5115-1
 Adaptive Fluidics™ Premium Pack



- Pack Includes:**
- | | | | |
|---|---|---|----------------------------|
| A | Cassette, tubing with pulse reflux and StableChamber® | E | Test chamber |
| B | Needle wrench | F | Bottle spike |
| C | 2.4mm yellow sleeves x2 | G | Tray cover drape |
| D | Plug | H | Screen drape |
| | | I | Remote control drape |
| | | J | MICS™ blue straight needle |



BL5115-2
 Adaptive Fluidics™ Standard Pack



- Pack Includes:**
- | | | | |
|---|------------------------------------|---|----------------------------|
| A | Cassette, tubing with pulse reflux | F | Bottle spike |
| B | Needle wrench | G | Tray cover drape |
| C | 2.4mm yellow sleeves x2 | H | Screen drape |
| D | Plug | I | Remote control drape |
| E | Test chamber | J | MICS™ blue straight needle |



BL5115-3
 Adaptive Fluidics™ Premium MICS™ Pack



- Pack Includes:**
- | | | | |
|---|---|---|----------------------------|
| A | Cassette, tubing with pulse reflux and StableChamber® | E | Test chamber |
| B | Needle wrench | F | Bottle spike |
| C | 1.8mm MICS™ sleeves x2 | G | Tray cover drape |
| D | Plug | H | Screen drape |
| | | I | Remote control drape |
| | | J | MICS™ blue straight needle |



BL5115-4
 Adaptive Fluidics™ Basic Pack



- Pack Includes:**
- | | | | |
|---|---------------------------|---|----------------------|
| A | Cassette, standard tubing | E | Test chamber |
| B | Needle wrench | F | Bottle spike |
| C | 1.8mm MICS™ sleeves x2 | G | Tray cover drape |
| D | Plug | H | Screen drape |
| | | I | Remote control drape |



*Only for Stellaris ELITE™ Anterior (2500 cpm)

Phaco Tips for MICS™ Surgery

RIS: Recommended Incision Size in Coaxial

OD: Outer Diameter

ID: Inner Diameter

Single Use

BL3318S MICS™ needle 1.8

STERILE   = 6



BL3318AS MICS™ needle 1.8 angled 30°

STERILE   = 6



RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm

BL3420S 2.0 MICS™ needle

STERILE   = 6

BL3420AS 2.0 MICS™ needle angled 30°

BL3118 Kit of 2 SU MICS™ Sleeves

STERILE   = 12




BL3124S Kit of 2 SU Sleeves 2.4mm
surgery



STERILE   = 12



Stellaris® System Coaxial I/A Handpieces

85794ST MICS™ I/A Handpiece, curved shaft with sleeve for 1.8mm incision STERILE   = 12



85795ST MICS™ I/A Handpiece, 45° angled tip with sleeve for 1.8mm incision STERILE   = 12



Stellaris® System CapsuleGuard® I/A Coaxial Handpieces

85914ST CapsuleGuard® I/A Handpiece MICS® Straight Tip - Incisions: sub 2.0mm STERILE   = 12



85915ST CapsuleGuard® I/A Handpiece MICS® 45° Tip - Incisions: sub 2.0mm STERILE   = 12



Stellaris® System Biaxial I/A Handpieces Aspiration

85901S19 Aspiration Handpiece 19G (Ø 1.07mm) STERILE   = 12

85901S20 Aspiration Handpiece 20G (Ø 0.89mm) STERILE   = 12

85780S Aspiration Handpiece 21G (Ø 0.8mm) STERILE   = 12

85790S Aspiration Handpiece 22G (Ø 0.7mm) STERILE   = 12

85792S Aspiration Handpiece 23G (Ø 0.6mm) STERILE   = 12



curved shaft roughened tip

Stellaris® System Biaxial I/A Handpieces Irrigation

85902ST19 1.8mm B-MICS 19G (Ø 1.07mm) Irrigation Handpiece STERILE   = 12

85787ST 21G Irrigation Handpiece (Ø 0.8mm) STERILE   = 12

85806ST 22G Irrigation Handpiece (Ø 0.7mm) STERILE   = 12

85807ST 23G Irrigation Handpiece (Ø 0.6mm) STERILE   = 12



curved shaft smooth tip

Stellaris® System Coaxial Zero Phaco Handpieces

SUP01ST Stellaris® System MICS™ Zero Phaco Handpiece 1.8mm STERILE   = 12



For aspirating softened nuclei during femtosecond cataract surgery. Designed with an open port aspiration tube with 30° bevel. Recommended for use in a 1.8 mm incision. Designed for use with the Stellaris® System.



SUP06ST Stellaris® System Zero Phaco Handpiece 2.2mm STERILE   = 12

For aspirating softened nuclei during femtosecond cataract surgery. Designed with an open port aspiration tube with 30° bevel. Recommended for use in a 2.2mm incision. Designed for use with the Stellaris® System.



SUP05 Zero Phaco Aspiration Handpiece 21G (0.8mm) STERILE   = 12

Bimanual aspiration handpiece, designed with an open port aspiration tube for aspirating softened nuclei during Femtosecond cataract surgery. 30° tapered tip. For use with 85787ST - Irrigation Handpiece 21G.





Product Reference	Product Name	Page
09.05.23PIN	Pick Forceps 23G	17
09.05.25PIN	Pick Forceps 25G	17
09.10.23PIN	Fine Tip Eckardt Forceps 23G	16
09.10.23PIN-M	Fine Tip Eckardt Forceps 23G, 37mm	17
09.10.25PIN	Fine Tip Eckardt Forceps 25G	16
09.10.25PIN-M	Fine Tip Eckardt Forceps 25G, 37mm	17
09.12.23PIN	End-Grasping Forceps 23G	16
09.12.25PIN	End-Grasping Forceps 25G	16
09.12PIN	End-Grasping Forceps 20G	16
10.24.23PIN	Morris Vertical Scissors 23G	17
10.24.25PIN	Morris Vertical Scissors 25G	17
10.24PIN	Morris Vertical Scissors 20G	17
11.20PIN	Foreign Body Forceps 19G	17
11.40PIN	BB Removal Forceps 18G	17
11.46.25PIN	Adeleberg Forceps 25G	17
12.03.25	Subretinal Cannula 25G, 41G Tip	18
12.06.25	Hubbard Visco Dissector 25G, 28G Tip	18
12.06.27	Hubbard Visco Dissector 27G	18
12.21PIN	Foreign Body Magnet 18G	17
20.04	Tano 20G	18
20.04.23	Tano 23G	18
20.04.25S (Stiff)	Tano 25G	18
20.04.27	Tano 27G	18
20.16	Retractable Tano 20G	18
20.16.23	Retractable Tano 23G	18
20.16.25	Retractable Tano 25G	18
23.76.23PIN	Horizontal Scissors Curved 23G	17
23.76.25PIN	Horizontal Scissors Curved 25G	17
23.76PIN	Horizontal Scissors Curved 20G	17
30.08.23	Soft Tip Cannula 23G	18
30.08.25	Soft Tip Cannula 25G	18
30.08.27	Soft Tip Cannula 27G	18
39.03.23PIN	Tano Asymmetrical Micro Forceps 23G	16
39.03.25PIN	Tano Asymmetrical Micro Forceps 25G	16
39.03.27PIN	Tano Asymmetrical Micro Forceps 27G	16
39.08.23PIN	Eckardt Forceps 23G	16
39.08.25PIN	Eckardt Forceps 25G	16
39.08.27PIN	Eckardt Forceps 27G	16
39.08PIN	Eckardt Forceps 20G	16
39.12.23PIN	Microserrated™ End-Grasping Forceps 23G	16
39.12.25PIN	Microserrated™ End-Grasping Forceps 25G	16
39.15.23PIN-M	Microserrated™ Tewari Omnicep Forceps 23G, 37mm	17
39.15.25PIN-M	Microserrated™ Tewari Omnicep Forceps 25G, 37mm	17
39.32.23PIN	Reddy End-Grasping forceps 23G	16
39.32.25PIN	Reddy End-Grasping forceps 25G	16
53.00.23	Bipolar Cautery 23G	18
53.00.25	Bipolar Cautery 25G	18
53.00.27	Bipolar Cautery 27G	18
55.25.23	Flexible Tapered Laser Probe 23G	11
55.25.25	Flexible Tapered Laser Probe 25G	11
55.26	Directional™ Laser Probe 20G	11
55.26.23	Directional™ Laser Probe 23G	11
55.26.25	Directional™ Laser Probe 25G	11
55.26.27	Directional™ Laser Probe 27G	11
55.47.23P	Illuminated Directional™ Laser Probe 23G	11
55.47.25P	Illuminated Directional™ Laser Probe 25G	11
55.70.23P	Flexible Tapered Illuminated Laser Probe 23G	11
55.70.25P	Flexible Tapered Illuminated Laser Probe 25G	11
65023S-12	Endgripping Forceps 23G	12

Product Reference	Product Name	Page
65025S-12	Endgripping Forceps 25G	12
85780S	Aspiration Handpiece 21G	22
85787ST	21G Irrigation Handpiece (Ø 0.8mm)	22
85790S	Aspiration Handpiece 22G	22
85792S	Aspiration Handpiece 23G	22
85794ST	MICS® I/A Handpiece, Curved Shaft	22
85795ST	MICS® I/A Handpiece, 45° Angled Tip	22
85806ST	22G Irrigation Handpiece (Ø 0.7mm)	22
85807ST	23G Irrigation Handpiece (Ø 0.6mm)	22
85901S19	Aspiration Handpiece 19G	22
85901S20	Aspiration Handpiece 20G	22
85902ST19	1.8mm B-MICS 19G (Ø 1.07mm) Irrigation Handpiece	22
85914ST	CapsuleGuard® I/A Handpiece MICS® Straight Tip	22
85915ST	CapsuleGuard® I/A Handpiece MICS® 45° Tip	22
BL14455	Stellaris ELITE™ without integrated laser	4
BL15455	Stellaris ELITE™ with integrated laser	4
BL2272	Laser Upgrade Kit	10
BL2296	Auxiliary Footpedal for LIO Use	11
BL2298	TruFocus LIO Premiere™ with LED illumination	11
BL2300HD	DMS (Digital Media System)	4
BL2394	Stellaris ELITE™ Integrated Foot Pedal	4
BL2394	Integrated Foot Pedal	10
BL3118	Kit of 2 SU MICS™ Sleeves	21
BL3124S	Kit of 2 SU Sleeves, 2.4 SURGERY	21
BL3170	Stellaris Phaco Handpiece	5
BL3234	Laser On/Off Keys (1 set)	10
BL3235	Laser Safety Goggles	10
BL3236	Zeiss Fixed ESF	10
BL3237	Leica Fixed ESF	10
BL3239	Laser Filter ESF Mount Topcon	10
BL3240	Zeiss 2 Position Filter ESF	10
BL3241	Leica 2 Position Filter ESF	10
BL3242	Adaptor for Fitting of 2 Position Filters	10
BL3244	Interlock Bypass Connector Kit	10
BL3245	Moeller Wedel Fixed ESF	10
BL3246	LIO Replacement Bulb	10
BL3247	Leica 4-way Filter	10
BL3270	Fragmentation Handpiece	5
BL3318AS	MICS® Needle 1.8 Angled 30° Single Use	21
BL3318S	MICS® Needle 1.8 Single Use	21
BL3320	Xenon Lamp Replacement	5
BL3321	Mercury Lamp Replacement	5
BL3394	Handle Extension	8
BL3420S	2.0 MICS® Needle	21
BL3420AS	2.0 MICS® Needle, Angled 30°	21
BL4363	Stellaris ELITE™ Zero Level Bottle Hanger - Dual Bottle	4
BL4364	Stellaris ELITE™/Stellaris Mayo Tray Adapter	5
BL4390	Foot Control Battery	4
BL4391	Foot Control Wall Charger	4
BL4393	Foot Control Charging Cradle	4
BL4394	Foot Control Power Cable	4
BL5115-1	Adaptive Fluidics™ Premium Pack	20
BL5115-2	Adaptive Fluidics™ Standard Pack	20
BL5115-3	Adaptive Fluidics™ MICS Premium Packs	20
BL5115-4	Adaptive Fluidics™ Basic Pack	20
BL5234	Adjustable Chandelier 23G	9
BL5235	Adjustable Chandelier 25G	9
BL5235S	Vivid Chandelier 25G	9
BL5237S	Vivid Chandelier 27G	9

Product Reference	Product Name	Page
BL5240	27G ESA Kit-Valved (kit of 3)	8
BL5241	27G Single Trocar-Valved	8
BL5275	20G Fragmentation Needle 15°	5
BL5277	23G Fragmentation Needle 30°	5
BL5280	Dual Infusion Kit	6
BL5281	23G ESA Kit-Valved (kit of 3)	8
BL5282	25G ESA Kit-Valved (kit of 3)	8
BL5283	23G Single Trocar-Valved	8
BL5284	23G Replacement Valves	8
BL5285	25G Single Trocar-Valved	8
BL5286	25G Replacement Valves	8
BL5290LC	Illuminated Laser Probe Curved 20G	11
BL5293LC	Illuminated Laser Probe Curved 23G	11
BL5295LC	Illuminated Laser Probe Curved 25G	11
BL5823M	Mid-Field 23G	9
BL5823D	Diffuse Wide Angle 23G	9
BL5825M	Mid-Field 25G	9
BL5825D	Diffuse Wide Angle 25G	9
BL5827D	Adjustable Chandelier 27G	9
BL5827M	Mid-Field 27G	9
BL7600	Universal Viscous Fluid Control Pack	19
BL7623	High Flow VFI Cannula 23G	19
BL7625	High Flow VFI Cannula 25G	19
CX9430	Reusable Bipolar Cord with Lemoso Connector	13
D31.08.23	23G Silicon Tip Active Backflush	18
D31.08.25	25G Silicon Tip Active Backflush	18
E7926	Bipolar Eraser 25G, Tapered, Lemoso Type Connector	13
E7928	Scleral Bipolar Eraser	13
E7929	Bipolar Eraser 20G	13
E7930	Bipolar Eraser 23G	13
E7932	Bipolar Eraser, Tapered, Lemoso type connector 27G	13
LP7519	LaseredgePlus® 19G MVR Blade	14
LP7520	LaseredgePlus® 20G MVR Blade	14
SE5423MVB/SE5423WVB	Bi-Blade® 25G Posterior Pack	7
SE5425MVB/SE5425WVB	Bi-Blade® 25G Posterior Pack	7
SE5427MVB	Bi-Blade® 27G Posterior Pack	7
SE5523MVB/SE5523WVB	Bi-Blade® 23G Combined Pack	6
SE5525MVB/SE5525WVB	Bi-Blade® 25G Combined Pack	6
SE5527MVB	Bi-Blade® 27G Combined Pack	6
SE5623B	Bi-Blade® 23G Vitrectomy cutter	8
SE5625B	Bi-Blade® 25G Vitrectomy cutter	8
SE5627B	Bi-Blade® 27G Vitrectomy cutter	8
SUI03-12	Asymmetric Peeling Forceps 23G	12
SUI05-12	Asymmetric Peeling Forceps 25G	12
SUI07-12	ILM Forceps 23G	12
SUI08-12	ILM Forceps 25G	12
SUI09-12	ILM Forceps 20G	12
SUI10-12	Serrated Forceps 20G	12
SUI11-12	Serrated Forceps 23G	12
SUI12-12	Serrated Forceps 25G	12
SUI13-12	SEndgripping Forceps 27G	12
SUI14-12	Serrated Forceps 27G	12
SUJ02	Passive Aspiration Handpiece. 1.5mm Soft Tip Cannula 25G	14
SUJ03	Passive Aspiration Handpiece. 1.5mm Soft Tip Cannula 23G	14
SUJ06	Passive Aspiration Handpiece. Blunt Needle 23G	14
SUJ07	Passive Aspiration Handpiece. Blunt Needle 25G	14
SUJ08	Backflush/Extrusion Handpiece. Blunt Tip 20G	14
SUJ09	Backflush/Extrusion Handpiece. Soft Tip Brush 20G	14
SUJ10	Backflush/Extrusion Handpiece. Soft Tip 20G	14

Product Reference	Product Name	Page
SUJ11	Backflush/Extrusion Handpiece. Blunt Tip 23G	14
SUJ12	Backflush/Extrusion Handpiece. Soft Tip brush 23G	14
SUJ13	Backflush/Extrusion Handpiece. Soft Tip 23G	14
SUJ14	Backflush/Extrusion Handpiece. Blunt Tip 25G	14
SUJ17	Backflush/Extrusion Handpiece. Blunt Tip 27G	14
SUK01	Membrane Pick 23G	13
SUK02	Membrane Pick 25G	13
SUP01ST	Stellaris® System MICS™ Zero Phaco Handpiece 1.8mm	22
SUP05	Zero Phaco Aspiration Handpiece 21G (0.8mm)	22
SUP06ST	Stellaris® System Zero Phaco Handpiece 2.2mm	22
VRL100	DK-Line® 5ml	19
VRL110	DK-Line® 7ml	19
VRL200	Okta-Line® 5ml	19
VRL500	Oxane® 1300, 10ml Syringe	19
VRL600	Oxane® 5700, 10ml Syringe	19
VRL700	Oxane® HD, 10ml Syringe	19

 @BauschSurgical

 Bausch + Lomb Surgical

 www.bauschsurgical.eu

Please contact your local representative for more information on Bausch + Lomb Retina products.

© 2021 Bausch & Lomb Incorporated.