

Stellaris Elite™

Vision Enhancement System



CATARACT



LASER



RETINA

BAUSCH + LOMB
See better. Live better.

EXPERIENCE THE DUAL-PORT DIFFERENCE

We are constantly seeking new ways to put innovation into your practice, and the Bi-Blade vitrectomy cutter is a perfect example of our dedication to providing precise and efficient surgical experiences.

STABILITY 15,000 CPM






















The Bi-Blade cutter's dual-edge blade effectively doubles the cut rate, which allows for faster procedures while potentially decreasing the risk of vitreoretinal traction ⁽¹⁾.

EFFICIENCY

Thanks to its unique dual port design that allows for cutting without port closure.

As the port remains open, the Bi-Blade provides consistent aspiration at all times during cutting. The Bi-Blade cutter's have been designed to create a consistent flow rate that operates independently from the cut rate.

1. Teixeira A. *et al.* Vitreoretinal traction created by conventional cutters during vitrectomy. *Ophthalmology* 2010;117:1387-1392

Stellaris Elite™ Combined Pack Bi-Blade®				
SE5525MVB / SE5525WVB 25G Combined Pack		 Sterile	 6	 Single Use
SE5527MVB 27G Combined Pack		 Sterile	 6	 Single Use
Stellaris Elite™ Posteriors Pack Bi-Blade®				
SE5425MVB / SE5425WVB 25G Posterior Pack		 Sterile	 6	 Single Use
SE5427MVB 27G Posterior Pack		 Sterile	 6	 Single Use
Stellaris Elite™ Bi-Blade® Vitrectomy Cutter				
SE5625B 25G Vitrectomy cutter Bi-Blade®		 Sterile	 6	 Single Use
SE5627B 27G Vitrectomy cutter Bi-Blade®		 Sterile	 6	 Single Use

EMEA_SU_SS_BB_18_001